

PROTECTOR

Hand Sanitizing Wipes

Form: Premoistened Towel

Color: Colorless

Odor: Mint

pH: N/A

Solubility (in water): Soluble*

(in mineral spirits): Insoluble

VOC Content (% by weight): N/D

Flash Point (ASTM D-7821): 70°F*

Specific Gravity: 0.88 g/cm³*

Density: 7.34 lbs/gal *

Storage Stability (at 70°F): 1 year

*Indicates test was performed on liquid component of towellette only

<u>This Product Contains:</u>	<u>CAS #</u>
Ethanol	64-17-5
Water.....	7732-18-5
Propylene Glycol	57-55-6
Linear Alcohol Ethoxylate.....	34398-01-1
Isopropyl Myristate.....	110-27-0

Great For Use In...

- **Medical Offices**
- **Hospitals**
- **Hotels**
- **Schools**
- **Dentist Offices**
- **Day Care Centers**
- **Restaurants**
- **Nursing Homes**
- **Correctional Facilities**
- **Emergency Vehicles**
- **Health Clubs**

DIRECTIONS FOR DISPENSING: Remove lid. Pull up the towel at the center of the roll by the corner and thread it through the “X” from the underside of the lid. Replace lid onto container and close hinged tab to prevent towels from drying out. When removing a towel, pull at a 90° angle until towel tears away at perforation.

DIRECTIONS FOR USE: Unfold a towel and thoroughly wipe down hands and skin for approximately 30 seconds or until towel and hands are dry. Repeat when necessary.

HMIS				NFPA	
		Severe	4	Extreme	
Health	0	Serious	3	High	Health..... 0
Flammability.....	3	Moderate	2	Moderate	Flammability..... 3
Reactivity.....	0	Slight	1	Slight	Reactivity..... 0
Personal Protection	None	Minimal	0	Insignificant	Special Precautions.....None